



Horizontal Laboratory Autoclave LHA-B13

www.labtron.com | info@Labtron.com

Overview

Horizontal Laboratory Autoclave LHA-B13 is a microprocessor controlled automatic autoclave. It has a capacity of 400 L and a maximum working temperature of 134 °C. Equipped with safety features like automatic shutoff preventing overload of current, auto protection against raised temperature and pressure making a reliable and user friendly device.

Features :

Microprocessor control with fully automatic cycles

Digital LCD display to monitor working status and parameters

Working pressure : 0.22 MPa

Working temperature : 134 °C

Automated drying function

Door lock feature : Door cannot be opened till the pressure inside the chamber doesn't fall to 0.027 MPa

Safe door lock system ensuring safety / protection during sterilization

Automatic power cut off with alarm with low water levels indication

Specifications :

Capacity	400 L
Sterilization Temperature	134 °C
Sterilization Pressure	0.22 MPa
Timer Range	0 - 60 minutes
Heat Average	$\leq \pm 1$ °C
Chamber Dimensions ($\Phi \times D$)	700 × 1100 mm
Product Dimension (L × W × H)	1430 × 880 × 1830 mm
Net Weight	420 kg
Gross Weight	530 kg
Power Supply	18 kW / AC 380 V 50 Hz

Applications :

Horizontal autoclave with cylindrical shape is suitable for medical, laboratory, microbiological, clinical use and to sterile the waste before discarding.



Quatro House, Lyon Way, Camberley, Surrey GU16 7ER United Kingdom

Email: info@Labtron.com | Website: www.labtron.com